

TISSUE Management

Trust.Worthy.Innovation.



DUX Dental Retraction Materials

Retraction Cord Dispenser

ShortCut™

(GingiBraid+ inside)



THE DENTAL
ADVISOR
++++ ½

Retraction Cord in a Bottle

GingiBraid+™

GingiKnit™

GingiGel™



THE DENTAL
ADVISOR
++++

Haemostatic Solution

Styptin™



Your Tissue Management Partner

ShortCut™


The innovative all-in-one retraction cord delivery system



ShortCut is available with GingiBraid+ retraction cord inside.

Dispensing cord with ShortCut is less cumbersome and definitely more efficient, convenient and hygienic than the cord in the bottle method.

ShortCut User Instructions



ShortCut allows hassle-free, scissor-free cord dispensing.

1. Simply turn the dial wheel clockwise to dispense the desired length of cord (each click 1 cm)
2. Squeeze the right side (cutter) for an effortless cut
3. Remove the cord

Posterior: 3 clicks / Anterior: 2 clicks.

Retraction cord delivery system



ShortCut (Braided Retraction Cord)

- All-in-one-system that fits comfortably in the hand
- Built-in cutter, no need for scissors
- Rotary mechanism dispenses each time the same amount of cord with a simple click (1cm)
- Durable, waterproof label
- Improved hygiene control

ShortCut is available non-impregnated or impregnated with 10% aluminum potassium sulfatet. **Sizes: 0 (thin), 1, 2 (thick)**

Retraction cords in a bottle



GingiBraid+ (Braided Retraction Cord)

Specially treated braid for faster absorption allowing immediate use and effectiveness.

- Will not split or collapse in the sulcus
- No memory effect
- Cords are easily identified by different color strands

GingiBraid+ is available non-impregnated or impregnated with 10% aluminum potassium sulfate. **Sizes: 0 (thin), 1, 2, 3 (thick)**



GingiKnit (Knitted Retraction Cord)

Constructed of many tiny, absorbent loops knitted in long, interlocking chains.

- Remains in-situ, even when touched by a bur or diamond
- Places easily without fraying or memory
- Cords are easily identified by different color strands

GingiKnit is available non-impregnated or impregnated with 10% aluminum potassium sulfate. **Sizes: 000 (very thin) 00, 0, 1, 2, 3 (thick)**




GingiGel (Gel-Coated Braided Retraction Cord)

This completely saturated braid is made of synthetic fiber and packed in a hemostatic gel.


- Provides immediate hemostasis and keeps the cord soft and pliable
- The cord does not dry out and is easily removed without subsequent bleeding

GingiGel is saturated and coated with 20% buffered aluminum chloride gel. **Sizes: 1 (thin), 2, 3 (thick)**


Effective Tissue Management - 5 steps to a successful retraction cord placement



Step 1:
The retraction cord is packed by positioning the instrument tip toward the starting point of the retraction cord to prevent inadvertent dislodging.



Step 2:
Insert gently below the free gingival margin to deflect the soft tissue away from the tooth.



Step 3:
Proper tissue management provides excellent access to the margin and an optimal final impression.



Step 4:
The retraction cord should be removed immediately before taking the impression.



Step 5:
Effective placement of retraction cord in either a single - or double-cord technique facilitates accurate impression making.

TISSUE Management



GingiBraid braided retraction cords are made of cotton fiber that allows fast absorption and will not split or collapse in the sulcus.



GingiKnit knitted retraction cords conform to sulcular anatomy with gentle, outward force and remain in situ even when touched by a bur.

Haemostatic solution



Styptin (Haemostatic Solution with Cranberry Flavour)

The aluminum chloride is a buffered 20% solution, which is quickly absorbed by the mucous tissue for fast haemostasis.

- Does not taste bitter thanks to the addition of cranberry flavour
- Blue coloring for improved visual control, without discoloring the gingiva

Styptin is supplied in a dispensing bottle (15ml) or a refill bottle (30ml).

Tissue Management order information

| Retraction Cords | Size | Length | Non-Impregnated* | Impregnated |
|--------------------------|------|--------|------------------|------------------|
| ShortCut (dispenser) | 0 | 183 cm | REF 13740 (0n) | REF 13730 (0a) |
| ShortCut (dispenser) | 1 | 183 cm | REF 13741 (1n) | REF 13731 (1a) |
| ShortCut (dispenser) | 2 | 183 cm | REF 13742 (2n) | REF 13732 (2a) |
| GingiBraid+ (bottle) | 0 | 183 cm | REF 13217 (0n) | REF 13210 (0a) |
| GingiBraid+ (bottle) | 1 | 183 cm | REF 13214 (1n) | REF 13211 (1a) |
| GingiBraid+ (bottle) | 2 | 183 cm | REF 13215 (2n) | REF 13212 (2a) |
| GingiBraid+ (bottle) | 3 | 183 cm | REF 13216 (3n) | REF 13213 (3a) |
| GingiKnit (bottle) | 000 | 183 cm | REF 13498 (000n) | REF 13508 (000a) |
| GingiKnit (bottle) | 00 | 183 cm | REF 13499 (00n) | REF 13509 (00a) |
| GingiKnit (bottle) | 0 | 183 cm | REF 13500 (0n) | REF 13510 (0a) |
| GingiKnit (bottle) | 1 | 183 cm | REF 13501 (1n) | REF 13511 (1a) |
| GingiKnit (bottle) | 2 | 183 cm | REF 13502 (2n) | REF 13512 (2a) |
| GingiKnit (bottle) | 3 | 183 cm | REF 13503 (3n) | REF 13513 (3a) |
| Haemostatic Solution | | Volume | Article number | |
| Styptin dispenser bottle | | 15 ml | REF 13157 | |
| Styptin refill bottle | | 30 ml | REF 13156 | |

*Impregnate with a gingival haemostatic solution such as **Styptin™** if desired.

Dux Dental

Zonnebaan 14 • NL-3542 EC Utrecht • The Netherlands

Tel. +31 30 241 0924 • Fax +31 30 241 0054

Email: info@dux-dental.com • www.duxdental.com