

Zone®

Trust.Worthy.Innovation.



TEMPORARY CEMENT

Zone hard setting crystalline cement combines superior, long-lasting retention with easy removal and clean up

Zone temporary cement is formulated completely free of eugenol. Zone cement will not soften, craze or discolour acrylic or polycarbonate crowns, is kind to tissue, nerve endings and pulp and will not decalcify dentin. Zone temporary cement does not leave a “rubbery” residue and can easily be removed.



Zone temporary cement does not leave a “rubbery” build up and is removed easily without time consuming and difficult cleaning!



Zone. Easy removal and clean up.



Zone® Temporary Cement

- Crystalline set cement combines excellent retention with easy removal
- Does not leave a “rubbery” residue on preparation or crown
- Eugenol free, non toxic
- Will not soften, crack or craze acrylics
- Does not discolor resin or polycarbonate crowns
- Does not decalcify the dentin
- Accepted by pulp and tissue
- Available in various dispensers

Dual Cartridge - Single Plunger Syringe

- Precise, even dispensation of base and catalyst every time
- Eliminates tubes and measuring
- One-handed dispensing

Ref 27043 (dual cartridge syringe 15 gr.)

Single Unit Dose

- Individual pouches for ultimate convenience
- Optimised asepsis
- Ideal for the patient who travels while provisionalized

Ref 27039 (25 pcs)

Zone Kit

- Kit of 2 tubes (base and catalyst)
- Including non-slip mixing pad

Ref 27040
(base 20 gr./catalyst 30 gr.)

Zone AutoMix Syringe

- Dispensed directly into any temporary
- Immediately seated
- Including 8 mixing tips
- Ultimate convenience: no mixing, no waiting, no hassle

Ref 27047 (syringe 4 gr.)

