


XENOLITE™ Technology

Trust.Worthy.Innovation.




LEAD-FREE X-RAY APRONS



Xenolite™ is a new hi-tech sheeting. Less weight, and just as effective as lead.



Easy to take on and off, with a comfortable thyroid-protection collar.



Both versions with Scotch-Guard™, durable non-slip corduroy backing.



Lead-Free Apron with collar (5 mm), Light Blue.



Panoramic Poncho Lead-Free (5 mm), Charcoal.

As light and as soft as a blanket, the Xenolite™ Lead-free X-Ray Aprons protect your patients from secondary radiation during intra-oral and panoramic radiography, without the inconvenience of lead. Covered with a new hi-tech alloy that ensures the same protection as lead, they weigh 30% less than standard lead aprons.



Xenolite. Combines comfort and protection!



• Lead-Free • Light Weight • High Tech Alloy (Xenolite™) • Panoramic Poncho

Provides protection for patients from secondary (scatter) radiation during intra-oral or panoramic radiography without the weight of lead. The Lead-Free X-Ray Apron and Panoramic Poncho feature a NEW high tech alloy sheeting that provides protection equivalent to a traditional Lead Apron, yet weighs 30% less! The Lead-Free X-Ray Aprons are much more comfortable for the patient, especially children, and are very easy for the technician to handle.

Features & Benefits:

- High tech alloy sheeting
- Easily cleaned micro fiber front
- Scotch-Guard™, durable non-slip corduroy backing
- Durable, stain-resistant binding
- Grommets permit hanging
- Protective multi-layered dental design
- Provides protection for patients from secondary (scatter) radiation during intraoral or panoramic radiography
- Inspected and certified by radiological lab
- Souple, drapes comfortably over patient; Thyroid collar wraps comfortably around neck
- Lead equivalency exceeds 0,5 mm
- Equal protection to standard lead / vinyl aprons
- Lightweight, easy handling
- Complies with known safe municipal landfill procedures
- Properly stored (draped, not folded), apron should last five to seven years



Product	Color	Ref
Lead free X-Ray Apron with collar (58 cm wide - 66 cm long)*	blue	31439
Lead free X-Ray Apron with collar, kids size (51 cm wide - 53 cm long)*	blue	31442
Lead free X-Ray Apron with collar (58 cm wide - 66 cm long)*	anthracite	31448
Lead free Panoramic Poncho (58 cm wide - 54 cm long = front side / 68 cm long = back side)**	blue	31440
Lead free Panoramic Poncho, kids size (51 cm wide - 36 cm long = front side / 51 cm long = back side)**	blue	31443
Lead free Panoramic Poncho (58 cm wide - 54 cm long = front side / 68 cm long = back side)**	anthracite	31449

* measured from neckline ** Shoulder measured from the centre of the neck